

Xevo G2-S/XS QTof LC-MS^E QC - UNIFI

Analytical Professional Services, Waters Corporation, Wilmslow, UK

BENEFITS

This document describes a SST (System suitability test). Routinely using this test will allow you to monitor system performance.

Critical parameters can be judged from comparison data like

- Solvent background level
- Stability of Retention time
- Mass Accuracy

WATERS SOLUTIONS

[Xevo G2-XS QTof](#)

[Unifi™](#)

EUanalytical_support@waters.com

KEYWORDS

Xevo G2-S/XS QTof, UNIFI, System Suitability Test, LCMS-QC

INTRODUCTION

This document is designed to help you test the complete LC-MS system prior to routine use. This test is designed so that stringent quality criteria (mass accuracy, RT stability, peak separation and solvent background) are assessed using a Waters certified standard solution (LCMS QCRM) in a 5 minute run. The method enables screening for 9 substances using the mass and retention time values as quality criteria. The method can also display and identify matched fragments, to further demonstrate the system's screening capabilities.

Waters recommends that you routinely perform this or a similar test to maintain the highest levels of confidence in your data and to assess the status of your instrument.

Checklist

- Check solvents and run ACQUITY LC system start-up to prime LC system and wash syringes
- Create a new analysis
- Go to initial conditions
- Run 3 Blanks and 6 LC-MS QC reference standard samples
- Review the report and ensure system meets performance criteria

Standards

Waters Part# [186006963](#) LCMS QC Reference Standard 1x 500 µL vial



LC conditions

| | |
|-----------------|--|
| Column: | 2.1 x 50mm ACQUITY BEH C18 1.7 µm |
| Column temp.: | 50 °C |
| Sample temp.: | 10 °C |
| Inj. volume: | 1 µL |
| Flow rate: | 500 µL/min. |
| Mobile phase A: | 0.1% formic acid in water 0.1% formic acid in acetonitrile |
| Mobile phase B: | 0.1% formic acid in acetonitrile |
| Wash | acetonitrile |
| Purge | 0.1% formic acid in water |
| Seal wash | 90:10 Water / methanol |
| Gradient | See Table |

MS conditions

| | |
|-------------------|--------------------------------------|
| System | Xevo G2-S/XS QTof |
| Ionisation mode | ESI ⁺ (ESI ⁻) |
| Acquisition mode | MS ^E , sensitivity mode |
| Acquisition range | <i>m/z</i> 50-1000 |
| Capillary voltage | 0.80 kV (< 2.0 kV)* |
| Collision energy | 20-40 eV ramping |
| Total run time | 5 min |

EXPERIMENTAL

LCMS QC Analysis parameters

Before you run the SST for the first time, import the analysis method, report template, sample list, custom fields, and library. Refer to the KCS article [WKB92598](#) "How to import the system suitability test materials for Vion IMS QTof?". Refer to *Table 1* for the Reference Compound details in order of elution.

Setup the LC system to the recommended LC conditions.

Use the following LC gradient for the analysis:

| Time (min) | Flow rate (ml/min) | %A | %B | Curve |
|------------|--------------------|----|----|---------|
| 0.00 | 0.500 | 95 | 5 | Initial |
| 0.50 | 0.500 | 95 | 5 | 6 |
| 1.50 | 0.500 | 75 | 25 | 6 |
| 2.00 | 0.500 | 75 | 25 | 6 |
| 3.50 | 0.500 | 30 | 70 | 6 |
| 3.90 | 0.500 | 30 | 70 | 6 |
| 4.00 | 0.500 | 95 | 5 | 6 |

Preparation of LCMS QC reference standard and blank sample

Prepare the lock mass and the calibration standard. Details for preparing the lock mass and the calibration solution can be found in KCS articles [WKB92649](#) and [WKB2458](#).

Run the Xevo G2-XS QTof LockSpray check (or LockSpray setup) and Calibration check (or Calibration). Ensure that fluidics bottle B contains sufficient lock mass solution.

Prime the LC system using the system start-up function.


Dilute 1:10 the LCMS QC reference standard in 95:5 water:acetonitrile + 0.1% formic acid. Prepare a blank sample with 95:5 water:acetonitrile + 0.1% formic acid.

| Component | Formula | Exact Mass [M+H] ⁺ | Exact Mass [M-H] ⁻ | Concentration µg/mL |
|--------------------|---|-------------------------------|-------------------------------|---------------------|
| Sulfaguanidine | C ₇ H ₁₀ N ₄ O ₂ S | 215.0597 | 213.0452 | 0.5 |
| Acetaminophen | C ₈ H ₉ NO ₂ | 152.0706 | 150.0561 | 1 |
| Caffeine | C ₈ H ₁₀ N ₄ O ₂ | 195.0877 | n/a | 0.15 |
| Val-Tyr-Val | C ₁₉ H ₂₉ N ₃ O ₅ | 380.2180 | 378.2034 | 0.25 |
| Leucine Enkephalin | C ₂₈ H ₃₇ N ₅ O ₇ | 556.2766 | 554.2620 | 0.25 |
| Sulfadimethoxine | C ₁₂ H ₁₄ N ₄ O ₄ S | 311.0809 | 309.0663 | 0.1 |
| Verapamil | C ₂₇ H ₃₈ N ₂ O ₄ | 455.2904 | n/a | 0.02 |
| Reserpine | C ₃₃ H ₄₀ N ₂ O ₉ | 609.2807 | 607.2661 | 0.06 |
| Terfenadine | C ₃₂ H ₄₁ NO ₂ | 472.3210 | n/a | 0.02 |

Table 1: LCMS QC mix content in order of elution


DATA ACQUISITION

From the My Work pane, click **Create analysis**, and then select **Create analysis by acquiring new data**.



Browse for the sample list (Xevo LC-MS QC SST) and analysis method (Xevo LC-HDMSE QC SST), specifying an analysis name and description.


From the newly created analysis tab, open the console pane by clicking the green tick on bottom right. Switch the MS to operate. From the summary menu, select "Go to initial conditions". The system should display "Idle", and the status "Ready". Wait for the Binary Solvent Manager delta pressure to be ≤ 20 PSI (≤ 1.37 bar), as shown in the figure below.



If necessary, amend the sample list vial position for the blank and the LCMS QC reference standard.

Click on the green start button to initiate the analysis.

Select "Run all samples" and make sure the options are ticked as depicted in the figure below



When the analysis is complete, click the **Report** tab and evaluate the results against the criteria defined in the "[Results Evaluation](#)" section. The report will use the following custom fields.

| Name | Description | Formula |
|-------------------------------|---|---|
| Average Peak Width in seconds | Average Chromatographic Peak Width in seconds of QC sample type | AVG(Chromatographic width,Sample type="QC")*60 |
| RMS ppm error | RMS ppm error of QC sample type | (SUM(Mass error^2,Sample type="QC")/COUNT(Mass error,Sample type="QC"))^0.5 |

RESULTS EVALUATION

After you have run the test, evaluate the results against the specifications in the table below:

| Criteria | Specitication |
|---|--|
| Mass accuracy | Better than 2 ppm RMS across all injections for all nine compounds |
| Retention time stability | Less than 3 seconds (0.05 minutes) RMS for the replicates |
| Peak width | Less than 3 seconds (0.05 minutes), measured as FWHM |
| BPI baseline level (solvent background) | Below 5% when evaluating the signal as XIC of the detected components in the "2: HD TOF MS" function (see figure 1 below). |

Item name: sst

Channel name: Identified Components




Figure 1: Typical Xevo G2-S/XS QToF LC-MS^E system suitability test chromatogram

If the mass accuracy acceptance criteria are not met, it is required to rerun the instrument MS setup. In case of poor chromatographic performance, prepare fresh mobile phases and inspect the LC flow path for potential leaks or blockages. Rerun the SST.

The SST can be run in negative ion mode by changing the ionisation polarity to negative ionisation as indicated in the [Experimental section](#). The targeted substances in the purpose tab must be re-imported to update the negative ion masses as indicated in [table 1](#). Note that caffeine, verapamil, and terfenadine are not detected in Negative Ion mode.